

Injection Trainer

NEW and IMPROVED! Multi-Layered soft tissue pad for subcutaneous and intramuscular injection.

Improved features include:

- More realistic skin that can withstand mulitple punctures without leaving puncture marks!
- Stronger muscle layer.
- More durable plastic base.
- Longer strap to fit most humans and manikins.

Skills

- Intramuscular injection
- Subcutaneous injection
- Professional to patient communication
- Managementof tissue

Features

- Multiple layers: dermis, fat & muscle
- Tissue is soft to the touch
- Liquid injected into the muscle layer can be squeezed out
- Trainer can be strapped onto arm, thigh, or stomach for role playing
- · Easy to use

Package supplied Injection Trainer

- 1 Skin and muscle pad
- 1 Strap on base
- 1 Injection needle

Injeciton Trainer - 3 Pack

• 3 Skin and muscle pads





Injection Trainer

- 3 Strap on bases
- 1 Injection needle